




Effectiveness of Carbonate Apatite Bone Graft Substitute by 1, 3, 4-wall bone defect model of Beagle dog

Yusuke Sakaguchi, Katsuyuki Yamanaka, Yutaka Shinozaki

GC Corporation



Introduction

Although an autograft is a gold standard of bone grafting, there is some limitation. Xenograft is also used but risk for infection potentially and we focused on chemically synthetic bone graft material. Focusing on bone composition, bone is a carbonate apatite containing carbonate ions, and it is known that carbonate ions play an important role in bone metabolism. As a result of extensive research, we developed a carbonate apatite bone graft substitute (Cytrans Granule, CO₃Ap, launched in Japan 2018) and a chemically synthetic resorbable membrane (Cytrans ElaShield, PBM, launched in Japan 2020).

In this study, we evaluated the effect of CO₃Ap and PBM on bone regeneration in vivo.

Composition of Bone




Figure 1. Composition of Bone

Inorganic composition of each material

Component	Contents (wt%)			
	Bone	HAp	β-TCP	Cytrans® Granules
PO ₄ ³⁻	46.6	56.8	61.2	48.0
Ca ²⁺	34.8	39.8	38.8	34.0
CO ₃ ²⁻	7.4	—	—	12.0
Na ⁺	0.9	—	—	4.15
Mg ²⁺	0.72	—	—	—
Cl ⁻	0.13	—	—	—
K ⁺	0.03	—	—	—
F ⁻	0.03	—	—	—



Figure 2. Cytrans® Granules
S(0.3-0.6 mm), M(0.6-1.0 mm), L(1.0-2.0 mm)



Figure 3. Cytrans® ElaShield
(Poly L-lactide- ε -caprolactone, PBM)

Method

The animal study protocol was approved by Institutional Animal Care and Use Committee of Hamri Co., Ltd., which has been approved by AAALAC International.


Study A:

After extraction of the bilateral mandibular anterior molars of beagle dog, it healed for 8 weeks. An 1-wall peri-implant bone defect (6 × 6 mm) was prepared and ϕ 3.0 mm dental implants were placed in distal side of the defect. The defect was filled with CO₃Ap(Cytrans Granules, GC). Sham group was set as a control. After a healing period of 12 weeks, μ CT and histological evaluation were performed. For histological evaluation, the newly formed bone area was measured by using microscopic images of undecalcified section (Villanueva Goldner stain) and Image J (NIH)




Study B:

A buccal dehiscence type 3-wall bone defect (10 × 4 × 5 mm) was prepared on the mandibular premolars. The defect was filled with CO₃Ap, and it was covered by PBM. A without-membrane group and sham group were set as a control. After 12 weeks of healing, then micro-CT and histological evaluation were performed. 3D reconstruction was performed using CT images. The volume of newly formed bone was then calculated using CT analysis software (VG Studio Max).



Study C:

Box-shaped defects (7 × 2 × 6 mm) were prepared, and ϕ 2.2 mm implants were placed in distal side of the defects. The defect was filled with CO₃Ap. Freeze-dried bone allograft (FDBA), deproteinized bovine bone mineral (DBBM), Hydroxy apatite (HAp) and β -Tricalcium phosphate (TCP) were set as controls. After healing period of 12 weeks, the samples along with the surrounding tissue were carefully harvested from the mandibular bone. Subsequently, Villanueva Goldner stain was applied to all sections and observed using a light microscope.



Results of Study A (1-wall)

In μ -CT analysis, the new bone volume was 139.0 mm³ at CO₃Ap, which was higher than 133.3 mm³ of the control. The new bone volume in the upper part of the defect was 12.0 mm³ at CO₃Ap and 7.7 mm³ at the control. CO₃Ap showed good regeneration especially in the upper part of the bone defect without loss of bone height. (Fig. 4-6)

As results of histological evaluation, new bone area was 46.3% in CO₃Ap group, which was higher than 38.5% in control group. (Fig. 7,8)




Fig 5. Bone volume using CT analysis




Fig 6. Bone volume in the upper part of defect using CT analysis




Fig 4. 3D CT image of mandibular bone (Left: Control, Right: CO3Ap, Yellow area: Defect area)




Fig 7. Histological image (Villanueva Goldner stain, Green; matured bone tissue, Red: osteoid)




Fig 8. Quantitative histomorphometry analysis of the new bone area

Results of Study B, C

Study B (3-wall) :

Histological evaluation revealed that the calcified area on CT observation was bone tissue in all cases. Newly formed bone volume ratio was 93.9% in the PBM group, 82.3% in the only CO₃Ap group, and 55.1% in the sham group.







	PBM+CO ₃ Ap	CO ₃ Ap	Sham
Inflammation	No	No	No
Micro CT image			
Histological Image (Villanueva Goldner Stain, Green: matured bone, Red: Osteoid)			

Figure 9. CT and Histological image




Figure 10. Newly formed bone volume by CT analysis

Study C (4-wall) :

As a result of histological evaluation, the descending order of bone augmented area was CO₃Ap, FDBA, DBBM, HAp and β -TCP. In FDBA, DBBM, HAp and β -TCP group, bone regeneration was insufficient at the upper part of the defect, however CO₃Ap group showed sufficient bone regeneration and maintained the bone height.




Figure 11. Histological image of each specimen (Osteoid : Magenta, Calcified bone : Yellow green)

Discussion

In this preliminary study, CO₃Ap showed the potential of being a superior and biocompatible bone graft substitute. CO₃Ap granules have been reported to be more susceptible to phagocytosis by osteoclasts than hydroxyapatite or β -tricalcium phosphate. It has also been reported that osteogenic differentiation of human bone marrow cells was accelerated with CO₃Ap compared to cell culture plates or HAp, suggesting that CO₃Ap facilitates osteogenic differentiation of bone marrow cells and osteoblasts¹⁾.

Furthermore, it promoted bone formation effectively when combined with PBM which was able to make sufficient space for the period required for bone tissue regeneration.

PBM is consist of two layers, porous layer and solid layer. Porous layer is for retaining bloods and solid layer is for inhibiting soft tissue from invading to bone defects. The two-layer structure provides a barrier property, strength, and space for bone regeneration.

1) M Fujioka-Kobayashi et al, J Biomed Mater Res, 2022, 110, 1524-1532

Conclusion

CO₃Ap is expected as clinically useful materials for bone graft treatment. Furthermore, it is thought that CO₃Ap can effectively regenerate bone by combining it with PBM.





Fig 12. Mechanism of Bone remodeling

Osteoconduction



Bone Substitution




Fig 13. CT image of Cytrans ElaShield