

Polyethyleneimine Stabilized Silver Nanoparticles as an Efficient and Selective Colorimetric Assay for Promethazine


Daim Asif Raja, Muhammad Raza Shah, Muhammad Imran Malik*

H.E.J. Research Institute of Chemistry, International Center for Chemical and Biological Sciences (ICCBS), University of Karachi, Karachi 75270, Pakistan


Email : daimasif.raja25@yahoo.com

Introduction

- **Promethazine** is FDA-approved first-generation antiemetic, antihistamine, and H1 receptor antagonist drug.
- Used for the treatment of vomiting and nausea caused by cancer chemotherapy, migraine episodes, narcotic therapy.
- It cause significant sedation, agitation, hallucinations, and possibly apparent life-threatening events or sudden infant death syndrome.



Promethazine



Branched PEI



Diff. color of AgNPs with their approx. sizes

- **Polyethylenimine (PEI)** is a synthetic polymer with high cationic charge density
- First used as gene delivery with great gene delivery potential.
- 25 Kda PEI considered as gold standard in gene Delivery.

- **Silver Nanoparticles (AgNPs)** are efficient in sensing purpose by signal enhancement.

- AgNPs possess good electrical and catalytic property


- Excellent biocompatibility.

- Lower cost

- Considering the present investigation, the colorimetric sensing was done for the Promethazine based on PEI (Mn= 25,000 g/mol) modified silver nanoparticles synthesized by chemical reduction method.

Synthesis and Characterization of PEI-AgNPs

- Ratio = PEI : AgNO₃ : 1 : 15




Continuous stirring




PEI-AgNPs

UV-visible Spectroscopy




(A) UV-visible spectrum of PEI-AgNPs

AFM



(B) 2-D AFM Image of PEI-AgNPs

Stability of PEI-AgNPs




(A) before and after exposure at 100°C (B) Spectra at different concentrations of NaCl


(C) pH variation

PEI-AgNPs and Drug Interaction

Spectroscopic Recognition of Promethazine




Calibration Curve




(A) UV-visible spectra by using various concentrations of PRO with PEI-AgNPs (B) Calibration curve for amount of PRO at 421 nm

- A simple and quick way to achieve stable PEI-AgNPs by introducing NaBH₄ as a reducing agent and PEI as a stabilizing agent.
- Colorimetric detection by PEI modified AgNPs is successfully established for Promethazine.
- Change in the morphology is confirm by using AFM, SEM, and DLS.
- Linear range of 0.005-200 μM, low LoD of 0.003 μM and LoQ of 0.008 μM was found.
- This sensor was successfully applied for the detection of PRO in environmental and biological sample.

FT-IR Analysis




AFM Analysis SEM Analysis Zeta size Zeta Potential




(A) (B) AFM images (II) SEM images (III) (IV) Zeta size and Zeta Potential of PRO with PEI-AgNPs

Reaction Mechanism




Schematic representation showing interaction between PEI-AgNPs with PRO

Interference Study



Real Sample Analysis



(A) Tap Water (B) Blood Plasma (C) Urine (D) Serum