


Rehabilitation of Maxilla with Implant Supported Fixed Detachable Prostheses & Mandible with RPD Using Locator R-Tx Attachment





Roya Zandparsa, DDS, MSc, DMD

Clinical Professor @ Tufts University School of Dental Medicine
CEO & Founder of Expert Dental Advisory

Rehabilitation of Maxilla with Implant Supported Fixed Detachable Prostheses



Pre-OP: Maxillary Interim Complete Denture




Four RP Internal Hex Implants Placed In Maxillary Arch

Rehabilitation of Maxilla with Implant Supported Fixed Detachable Prostheses



Maxillary Interim Complete Denture over Titanium Healing Caps



**Integrated implants
with Titanium Healing
Abutments**




**Maxillary PMMA Interim Hybrid
Prostheses with Titanium
Temporary Copings**



**Monolithic Zirconia Hybrid
Prostheses with Definitive Straight
& Angled Tapered Abutments**

Definitive Monolithic Zirconia Maxillary Hybrid Prostheses



Implant-Supported Mandibular RPD



**Mandibular Arch with
Implant Site 29-30**




R-Tx Attachment



**Denture Attachment
Housing Cap with Black
Processing insert**



**Existing Mandibular RPD
With attachment**



**LOCATOR R-Tx Abutment 4.5 mm
on Internal Hex Implant**



**Locator R-Tx
Removable Attachment System**