

Active and Passive Mineralization of Bio-Gide® Membranes in Rat Calvaria Defects

Karol Ali Apaza Alccayhuaman DDS, PhD^{1,2}, Ulrike Kuchler DDS, PhD³ and Reinhard Gruber Dr. nat. Techn., PhD


¹ Department of Oral Biology, University Clinic of Dentistry, Medical University of Vienna; ² Karl Donath Laboratory for Hard Tissue and Biomaterial Research, University Clinic of Dentistry, Medical University of Vienna, ³ Department of Oral Surgery, University Clinic of Dentistry, Medical University of Vienna.




Background



In clinical cases with severe bone loss, we opt for GBR as our preferred clinical method.



For several years, collagen membrane (CM) has been used only as a mechanical barrier in combination with bone grafts¹.



However, recent studies showed CM alone exhibits osteoconductive properties and can lead to spontaneous mineralization^{2,3}.

1. Mizraji, G.; Davidson, A.; Gursoy, M.; Gursoy, U.K.; Shapira, L.; Wilensky, A. Membrane barriers for guided bone regeneration: An overview of available biomaterials. *Periodontol* 2000 **2023**,.

2. Feher, B.; Apaza Alccayhuaman, K.A.; Strauss, F.J.; Lee, J.S.; Tangl, S.; Kuchler, U.; Gruber, R. Osteoconductive properties of upside-down bilayer collagen membranes in rat calvarial defects. *Int J Implant Dent* 2021, **7**, 50,

3. Alccayhuaman, K.A.A.; Tangl, S.; Blouin, S.; Hartmann, M.A.; Heimel, P.; Kuchler, U.; Lee, J.S.; Gruber, R. Osteoconductive Properties of a Volume-Stable Collagen Matrix in Rat Calvaria Defects: A Pilot Study. *Biomedicines* 2021, **9**,




Hypothesis and Aim

Hypothesis: The use of minced CM will compromise the integrity of the membrane, thus disrupting its osteoconductive properties and spontaneous active mineralization.

Aim: We aimed to shift active mineralization in intact collagen membranes towards passive mineralization.

Material and Methods



Results





Micro- CT Analysis

- 3D images depicted mineralized bone volume and bone islands in the minced group, which looked like real bone due to strong density.
- Bone Volume (BV), trabecular thickness (Tb.Th), and coverage in the minced group almost reached values similar to those in the intact group.

Intact CM




Minced CM



Histology


Intact CM

The image shows an intact CM(A), with new bone growth visible within it (B) and towards the center (C). This new bone contains collagen fibers from the membrane named hybrid bone (D).



Minced CM

The image shows a minced CM that has mineralized areas (A). There are clusters of fibers (B) resembling bone tissue, but no bone cells or osteoid matrix (C) were found. The mineralized collagen fibrils interconnect the membrane fragments (D).



Conclusion

Collagen membranes (CM) allow biomineralization, independent of bone cells and osteoid matrix. This is promising because we can exploit the intrinsic properties of CM to accelerate bone repair and regeneration in the clinical setting.