




Novel plasma biomarkers track cognitive performance using digital testing

Monday-0364

Toniolo S^{1,2}, Zhao S³, Broulidakis MJ³, Gendarini C^{1,4}, Tan HK J^{5,6}, Sinha D^{5,6}, Tahira N⁵, Fower A¹, van Duijn C⁷, Thompson S², Manohar S^{1,2,3}, Koychev I⁴, Husain M^{1,2,3}

- Acetylcholinesterase (ACHE) levels were the best predictors of performance at digital testing.
- Novel markers such as PRDX6, IL-6, and VEGFD had a mild protective effect on cognition.
- Measures of long-term memory were the best predictors of ACHE and pTau217 levels.

INTRODUCTION	METHODS	RESULTS	RESULTS	CONCLUSIONS
Plasma biomarkers can identify the presence of Alzheimer's disease (AD) in clinical populations in a rapid and cost-effective way. However, measuring cognition remains essential. Digital tests can detect early signs of cognitive impairment at scale. At present, studies looking at the relationship between plasma biomarkers and performance at digital testing have looked only at a small array of biomarkers (e.g. pTau217, GFAP, NEFL). 391 participants (subjective cognitive decline (31), mild cognitive impairment (28), AD dementia (88), Lewy body disease (18), Corticobasal syndrome (9), Frontotemporal dementia (17) and age-matched controls (201)) were recruited. They were tested on a visual short-term memory task (Oxford Memory Test (OMT)), a fully remote online cognitive assessment tool part of Oxford Cognitive Testing Portal (OCTAL: https://octalportal.com), (Figure 1).	Three key metrics were computed: Accuracy (proportion of correctly identified items), Reaction Times and Localisation Error (distance between the original item and the response location).  Figure 1 Oxford Memory Test (OMT) task design. From the OCTAL platform the Oxford Memory Test was used in the full sample. Participants were presented with either 1 or 2 objects randomly distributed in the screen, after a second delay two objects appeared at the corners of the screen, one of which had appeared in the memory array whereas the other one was a distractor. Firstly, they needed to identify the object they had seen previously and then drag it back to its original location. Metrics extracted from the task included: Accuracy (proportion of correctly identified items), Reaction Times (time to drag the item to its remembered location) and Localisation Error (distance between the original item and participants' response location).	Plasma Acetylcholinesterase (ACHE) levels were the best predictors of performance at digital testing across the three short-term memory metrics studied (Figure 2).  Figure 2 Performance of key plasma biomarkers in predicting cognitive impairment at digital testing. Plasma Acetylcholinesterase (ACHE) levels were the best predictors of performance at digital testing across the three different metrics studied (Accuracy, Reaction Times, Localisation Error). For visualization purposes accuracy represents the inverse of the proportion between correctly identified/incorrectly identified items. Positive associations are shown in red, negative in blue, non-statistically significant associations in black.	ACHE levels were significantly predicted by all cognitive domains (panel a), were increased in AD dementia patients compared to controls, and particularly in patients who were on Acetylcholinesterase inhibitors, unlike pTau217 (panel b).  Figure 3 Comparison between ACHE and pTau217 levels' prediction by cognitive domains, and effects of medications. Panel a: Long-term memory measures were the best predictors of ACHE and pTau217 levels. ACHE levels were predicted by all cognitive domains, unlike pTau217. For visualization purposes all the metrics were set as higher values corresponding to worse performance. Positive associations are shown in red, non-statistically significant associations in black. A measure of the size of the association is available at the bottom of the plot. LTM = long-term memory, calculated from the delayed recall measure of the Object-in-Box Task; STM = short-term memory, corresponding to Accuracy measure of the OMT task; Vocabulary = vocabulary abilities measured by a digit span version of the Rey-Osterich Complex Figure Test; Processing Speed = reaction times of the Digit Symbol Substitution Test. Panel b: ACHE levels were increased in AD patients compared to controls, and particularly in patients who were on Acetylcholinesterase inhibitors (AChEi), unlike pTau217.	Studying a wider array of plasma biomarkers is beneficial for predicting cognitive performance. ACHE levels seem promising in measuring cognitive outcomes at digital testing. Further studies should unravel the best combination of plasma and digital biomarkers in patient stratification. Sofia Toniolo Senior Research Fellow and Honorary Consultant Neurologist Email: sofia.toniolo@ndcn.ox.ac.uk